




Biopsy Curettes

Malleable



04-6265
ASCH
Uterine
Secretion scoop
200 mm

Malleable



04-6266
Uterine
Secretion scoop
270 mm



04-6267
NOVAK-SCHOECKAERT
Endometrial biopsy
curette
240 mm



04-6268
RANDALL
Endometrial
biopsy curette
240 mm




Biopsy Forceps






Biopsy Forceps



- 2 Different shapes of the jaw are available
- a) Gothic shape
2.0 resp 2.4mm wide (causes less bleeding)
- b) Rectangular shape, 2.4mm wide



Biopsy Forceps






Uterine Scoops


Table of Recamier		Sharp		Blunt
Fig	Width	Rigid	Malleable	Rigid
00	4.5mm	04-6294		
0	7.5mm	04-6295		04-6306
1	8.5mm	04-6296	04-6317	04-6307
2	10.5mm	04-6297	04-6318	04-6308
3	12.0mm	04-6298	04-6319	04-6309
4	14.0mm	04-6299	04-6320	04-6310
5	15.0mm	04-6300	04-6321	04-6311
6	16.5mm	04-6301	04-6322	04-6312
7	18.0mm	04-6302	04-6323	04-6313
8	19.5mm	04-6303	04-6324	04-6314
9	21.5mm	04-6304		04-6315
10	23.5mm	04-6305		04-6316
04-6294- 04-6305		Sharp - Rigid		
04-6317- 04-6324		Sharp - Malleable		
04-6306- 04-6316		Blunt - Rigid		



04-6292
SIMON
max.240 mm



04-6293
SCHRODER
max.310 mm



04-6294
RECAMIER
max.300 mm






Fig	Width	Sharp		Blunt
		Rigid	Malleable	Rigid
1	7.5mm	04-6325	04-6331	04-6337
2	8.5mm	04-6326	04-6332	04-6338
3	10.5mm	04-6327	04-6333	04-6339
4	12.0mm	04-6328	04-6334	04-6340
5	14.0mm	04-6329	04-6335	04-6341
6	15.0mm	04-6330	04-6336	04-6342
04-6325 - 04-6330		Sharp - Rigid		
04-6331 - 04-6336		Sharp - Malleable		
04-6337 - 04-6342		Blunt - Rigid		

Fig	Width	Blunt - Malleable
0	16.5mm	04-6343
1	18.5mm	04-6344
2	19.5mm	04-6345
3	21.5mm	04-6346
4	23.5mm	04-6347
5	26.0mm	04-6348
6	28.5mm	04-6349

BUMM





Placenta-and Ovum Curettes



04-6351
COLLIN
300 mm



04-6352 - 04-6353
RECAMIER
blunt-rigid
max. 360 mm



Fig:1



Fig:2



04-6354 - 04-6358
BUMME
max. 310 mm
04-6354 - 04-6355 sharp-malleable
04-6356 - 04-6358 blunt-malleable



Fig:1



Fig:2




Fig:3

Fig	Width	Sharp Malleable	Blunt Malleable
1	23.5mm	04-6354	04-6356
2	26.0mm		04-6357
3	28.5mm	04-6355	04-6358



Vacuum Curettes

04-6359 - 04-6364

Suction curettes, stainless steel, fitting hose connection.

Working tips are smoothly rounded off, therefore no danger of perforation.



(suitable for suction biopsy)

04-6359: diam 4 mm

04-6360: diam 6 mm

04-6361: diam 8 mm

04-6362: diam 10 mm

04-6363: diam 12 mm

04-6364: diam 14 mm

04-6365

Hose connection inside diam. 11 mm, curved, stainless steel.



OLIVE
Ø 11.15mm

04-6365

04-6366

The new connection between suction curette and hose connection, (04-6365) Guarantees quick changing and free turning of the suction curettes on the hose connection.

04-6366

Transparent suction hose 1.8 m long, inside diam 11 mm, outside diam 16.5 mm can be sterilized in an autoclave at 1 bar overpressure (120 °C)



04-6367

Strainer for tissues, plastic, for single use only, to be used for histopathological examination.



Gynaecological Instruments



04-6368 diam 8 mm
04-6369 diam 10 mm
BOZEMANN-FRITSCH
Uterine Flushing Tube
350 mm



04-6371



04-6372

04-6370

SCHULZE

Salpingograph, for utero-salpingography
in wooden case, consisting of:

- 04-6371** Manometer only
- 04-6372** Connection only, without manometer
- 04-6373** Uterine cannula with tape without cone
- 04-6374** Cone only, small
- 04-6375** Cone only, medium
- 04-6376** Cone only, large
- 04-6377** Record syringe, 10CC, luer-lock mount.



04-6373



04-6374



04-6375



04-6376

04-6377



04-6378

SCHULZE

Salpingograph, equipment like(04-6370)
but without manometer.

