



Vacuum Curettes for Aspiration Lipectomy, Suction Instruments



07-9043 - 07-9048
BIEMER

07-9043 ∅ 4 mm
07-9044 ∅ 6 mm
07-9045 ∅ 8 mm
07-9046 ∅ 10 mm
07-9047 ∅ 12 mm
07-9048 ∅ 14 mm

Curettes for aspiration
lipectomy, graduated,
slightly curved



07-9049 - 07-9050
BIEMER

07-9049 ∅ 4 mm
07-9050 ∅ 6 mm
Curettes for aspiration
lipectomy, with three
apertures, slightly curved
malleable



07-9051 - 07-9053
BIEMER

07-9051 ∅ 8 mm
07-9052 ∅ 10 mm
07-9053 ∅ 12 mm
Curettes for aspiration
lipectomy, with three
apertures, slightly curved
malleable



07-9054

For aspirating
fluid from the excretory
ducts of large glands
and for keeping the
operating area clear
170 mm



07-9055

Trocar
For puncturing
150 mm



07-9056

For clearing
large quantities
of fluid
190 mm



Suction Instruments

07-9057 150mm

Suction handle, with vacuum regulator and finger cut-off, brass, mirror-finish chromium-plated



olive
Ø 9 mm

07-9058 170mm

Suction handle, with hand regulation valve and quick-shutter, brass, mirror-finish chromium-plated



Both suction handles have bayonet fitting olive and can be used for the suction cannulae. Ø 9 mm



07-9059
Trocars
For puncturing
170 mm



07-9060
Hollow needle,
sharp
150 mm



07-9061
double curved
For general
surgery
230 mm



07-9062
For tonsillectomy,
oral and laryngeal
operations
190 mm




07-9063
For tonsillectomy,
oral and laryngeal
operations
220 mm

★ With bayonet fitting to be used with suction handles (07-9057 - 07-9058)



Suction Instruments






Suction Instruments






Suction Instruments





Suction Instruments





Combined Suction and Irrigation Cannulae, Coagulation Suction Cannulae



07-9123
FINSTERER
320 mm

olive $\varnothing$ 9-13 mm $\frac{1}{2}$

$\varnothing$ 2 mm


$\varnothing$ 4,5 mm

07-9124
YASARGIL
LUER-LOCK cone
200 mm



07-9125
with finger cut-off and stylet,

olive 4,5 mm $\frac{1}{2}$



07-9126 - 07-9128
HOUSE

	suction tube diam.	irrigation tube Diam.
07-9126	0,9 mm	0,7 mm
07-9127	1,2 mm	0,9 mm
07-9128	1,5 mm	1,2 mm

with finger cut-off and stylet,
with LUER hub

$\varnothing$ 4 mm
olive
olive $\varnothing$ 9-6 mm $\frac{1}{2}$



07-9129 3,5 mm
with LUER-LOCK cone

olive $\varnothing$ 9 mm $\frac{1}{2}$



07-9130 $\varnothing$ 5.0 mm
LANDOLT
with finger cut-off and
stylet, can be sterilized
in an autoclave up to 3 bar (143 °C)

Connection for
suction hose

$\varnothing$ 7,5-12 mm
olive

4 mm diam., Untarnished
contact to be touched by
HF-electrodes